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
**Table 1. Choose the Quant-iT™ Assay Kit that meets your needs.**

Product	Sample selectivity	Applications
Quant-iT™ dsDNA HS Assay Kit, *0.2–100 ng*	dsDNA samples ranging from 10 pg/μl to 100 ng/μl	Genomic DNA Viral dsDNA Rare DNA PCR products
Quant-iT™ dsDNA BR Assay Kit, *2–1000 ng*	dsDNA samples ranging from 100 pg/μl to 1000 ng/μl	Plasmid preps Subcloning fragments DNA for transfections
Quant-iT™ RNA Assay Kit, *5–100 ng*	RNA samples ranging from 250 pg/μl to 100 ng/μl	Microarray samples Real-time PCR samples
Quant-iT™ Protein Assay Kit, *0.25–5 μg*	Protein samples ranging from 12.5 μg/ml to 5 mg/ml	PAGE gels Western blots Activity assays

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
**Figure 2. Performance of the Quant-iT™ RNA Assay Kit.**


The Quant-iT™ RNA Assay Kit has a linear detection range of 5–100 ng and is selective for RNA, even in the presence of an equal mass of DNA. The x-axis gives the mass of nucleic acid at a given point when DNA or RNA is assayed alone. In the 1:1 mixture, the total mass of nucleic acid at a given point, therefore, is double what is stated on the x-axis.

**Figure 3. Performance of the Quant-iT™ Protein Assay Kit.**

The Quant-iT™ Protein Assay Kit has a detection range of 0.25–5 μg. This graph shows the signal generated from nine different proteins over the stated detection range. The protein-to-protein variability at the 3 μg point for all nine proteins assayed is 12%.

**Figure 1. Performance of the Quant-iT™ High-Sensitivity DNA Assay Kit.**

The Quant-iT™ High-Sensitivity DNA Assay Kit has a linear detection range of 0.2–100 ng and is selective for dsDNA, even in the presence of an equal mass of RNA. The x-axis gives the mass of nucleic acid at a given point when DNA or RNA is assayed alone. In the 1:1 mixture, the total mass of nucleic acid at a given point, therefore, is double what is stated on the x-axis.

Figure 4. The Qubit™ fluorometer.

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Para mayor información dirigirse a Invitrogen Argentina S.A.